

The Development of Low Flow Pump Tubes & their Potential as a Tool for Simplifying the Validation of Axial Diffusive Tubes

The Development Of Low-Flow Pump Tubes And Their Potential As A Tool For Simplifying The Validation Of Axial Diffusive Tubes



Nicola Watson
Markes International

MARKES
International

www.markes.com

Overview

- Introduction to diffusive sampling including principles of diffusion
- Methods to validate uptake rates
- Diffusion locking technology and its use to calculate/confirm uptake rates for sorbent/analyte combinations.
- Parallel study results



MARKES
International

www.markes.com


Sorbent Tubes - Diffusive Sampling

- A simple and cost effective method of collecting the large number of samples required in many air monitoring programmes
- Vapors migrate across the air gap at a constant "uptake rate"
- Uptake rates for common compounds and sorbents available in the literature
- Diffusive sampling is a slow process, typically sample for days
- Only single bed tubes can be used



Brass Cap

Diffusion Cap (fitted to grooved end of tube)



MARKES
International

www.markes.com

Human Exposure Monitoring

1,3-Butadiene and benzene are criteria air pollutants cited in many regulations


Diffusive sampling facilitates personal exposure monitoring

Butadiene

- Pumped sampling of butadiene onto Carbopack X or Carbograph 5TD
- Diffusive sampling on Carbopack X:
 - 8 hour: 1.64 pg/ppb/min
 - 1 week: 1.19 pg/ppb/min
 - (2 week: 1.02 pg/ppb/min)

Benzene

- Validated for pumped and diffusive sampling at ambient and workplace levels. Optimum sorbent Carbograph 1TD.



MARKES
International

www.markes.com

Principles of Diffusive Sampling


- Diffusive Monitors must be capable of maintaining the following conditions during sampling
 - Ambient concentration of the analyte at the surface of the monitor
 - Zero concentration of the analyte at the surface of the sorbent
 - A linear concentration gradient between the two
- When these conditions are met, Fick's First Law of Diffusion applies

$$J = -D \frac{\partial \phi}{\partial x}$$

MARKES
International

www.markes.com

Axial Diffusive Sampling



- Using an industry standard tube (typical diffusion path length of ~1.5 cm) atmospheric concentration of a compound is determined by using the following equation.

$$\text{ppm} = \frac{\text{Mass of sample on tube (ng)}}{\text{Up (ng/ppm/min)} \times \text{T}_s \text{ (min)}}$$


MARKES
International

www.markes.com

The Development of Low Flow Pump Tubes & their Potential as a Tool for Simplifying the Validation of Axial Diffusive Tubes

Axial Diffusive Sampling


- Many actual uptake rates for a variety of compounds are now published in the literature.
- However, if an uptake rate is not available then there are several options,
 - Calculate ideal uptake rate using Fick's 1st law.
 - Predict the uptake rate by combining simple short experiments and modelling studies.
 - Determine the uptake rate experimentally.



MARKES International

www.markes.com

Specialist Safe-Lok Sorbent Tubes



- Why use low flow pumped tubes ?
 - Validation of uptake rates using an active flow similar to the natural rate of diffusion, normal axial tubes have errors associated with diffusion
 - Parallel sampling of air using diffusive tubes and low-flow pumped tubes can be a way to validate uptake rates

MARKES International

www.markes.com

Sorbent Tubes – Specialist tubes

- Safelok tubes facilitate pumped sampling at low flow rates (< 1 ml/min)
- Reduces risk of contamination
- Prevents necking caused by over tightening
- Safer to handle toxic compounds
- Same mass of sorbent and same external dimensions as standard tubes



MARKES International

www.markes.com

Experimental

- Tube types
 - Safelok tubes: Carbograph 1TD/Carbopack-X
 - Diffusive tubes: Carbopack-X
 - Safelok tubes pumped at 2 ml/min for 14 days
 - Diffusive tubes exposed for 14 days
- Analysis conditions
 - Tube desorption temp: 400°C for 15 min
 - Flow path temp: 200°C
 - Cold trap low: -30°C
 - Cold trap temp high: 300°C for 2 min
 - Trap Split: 18 ml/min
 - Cold trap packing: TENAX GR
 - Column pressure: 12.2 psi
 - Column used OV-1701
 - GC temp program 50°C for 6 mins then 10°C/min to 90°C then 10°C/min to 200°C for 10mins (31 mins total)



MARKES International

www.markes.com

RESULTS FROM THE ANALYSIS OF SAFELOK TUBES

Diffusive Tubes							
Tube	Duration	Benzene	Toluene	Ethylbenzene	m/p-Xylene	o-Xylene	Total Hydrocarbons
	min	ppb v/v	ppb v/v	ppb v/v	ppb v/v	ppb v/v	ppb v/v
CX-90	BLK	18991	<0.1	<0.3	<0.3	<0.3	<5
CX-537		18991	0.5	1.1	<0.3	0.5	10
CX-587		18991	0.4	1.0	<0.3	0.5	8
CX-3750		18991	0.4	0.9	<0.3	0.5	10
CX-502		18991	0.4	1.0	<0.3	0.5	9

Pumped Tubes							
Tube	Volume	Benzene	Toluene	Ethylbenzene	m/p-Xylene	o-Xylene	Total Hydrocarbons
	ml	ppb v/v	ppb v/v	ppb v/v	ppb v/v	ppb v/v	ppb v/v
M-029485	BLK		<0.1	<0.3	<0.3	<0.3	<5
M-029480	P37	34140	0.4	1.13	<0.3	0.5	7
M-029489	PL7	33094	0.4	1.03	<0.3	0.5	7
M-02481	P46	30591	0.4	1.07	<0.3	0.5	7
CX-140	P37 BU	34140	<0.1	<0.3	<0.3	<0.3	<5
CX-102	PL7 BU	33094	<0.1	<0.3	<0.3	<0.3	<5
CX-3253	P46 BU	30591	<0.1	<0.3	<0.3	<0.3	<5

MARKES International

www.markes.com

Calculating uptake rates

- For example Benzene has a concentration of 0.4 ppb on the pumped tube which also relates to 0.4 ppb on the diffusive tube.
- Given this, we can now calculate the uptake rate by




$$\text{Uptake rate} = \frac{\text{Mass (pg)}}{\text{Concentration (ppb)} \times \text{Time (min)}}$$

Which gives a value of 1.99 pg.ppb⁻¹.min⁻¹ for benzene on carbopack-x

MARKES International

www.markes.com

The Development of Low Flow Pump Tubes & their Potential as a Tool for Simplifying the Validation of Axial Diffusive Tubes



Summary

- Demonstrated how low flow pumped sampling can be used as a confirmation of uptake rates for variety of sorbent analyte combinations.
- Diffusion locking technology facilitates low flow sampling by minimising the error related to simultaneous diffusive ingress.
- Acknowledgements
 - Malcolm Henderson and Bob Lipscombe – NPL, Teddington

MARKES International www.markes.com